

CORELINE

PRODUCT CATALOG



CORELINE™

 **Hanger**
SUPPLY CHAIN



APOLLO AIR WALKING BOOT

L-Code: L4361 / L4360

- Built in air pump
- Rocker bottom provides natural gait
- Circumferential design for maximum support
- Removable toe guard

Size	Men's Shoe Size	Women's Shoe Size	Tall Product No.	Short Product No.
XS	<4	<5.5	CLTLXS	—
SM	4.5 - 7	6 - 8	CLTLSM	CLSHSM
MD	7.5 - 10	8.5 - 11.5	CLTMD	CLSHMD
LG	10 - 13	11 - 15	CLTLG	CLSHLG
XL	13+	15+	CLTLXL	CLSHXL



CORELINE NON-PNEUMATIC FULL SHELL TALL

L-Code: L4387

- Rocker bottom provides natural gait
- Circumferential design for maximum support
- Low profile design

Size	Men's Shoe Size	Women's Shoe Size	Product No.
XS	<4	<5.5	CLFSTXS
SM	4.5 - 7	6 - 8	CLFSTSM
MD	7.5 - 10	8.5 - 11.5	CLFSTMD
LG	10 - 13	11 - 15	CLFSTLG
XL	12+	15+	CLFSTXL



CORELINE NON-PNEUMATIC EU WALKER

L-Code: L4387 / L4386

- Rocker bottom base designed to enable a natural and stable gait
- Shock absorbing sole reduces impact on the heel while walking to prevent irritation of injuries

Size	Men's Shoe Size	Women's Shoe Size	Tall Product No.	Short Product No.
XS	<4	<5.5	CLTWXS	CLSWXS
SM	4.5 - 7	6 - 8	CLTWSM	CLSWSM
MD	7.5 - 10	8.5 - 11.5	CLTWMD	CLSWMD
LG	10 - 13	11 - 15	CLTWLG	CLSWLG
XL	13+	15+	CLTWXL	CLSWXL



CORELINE PNEUMATIC EU WALKER

L-Code: L4361 / L4360

- Features an air bladder system that inflates around the ankle to provide customized compression

Size	Men's Shoe Size	Women's Shoe Size	Tall Product No.	Short Product No.
XS	<4	<5.5	CLTWAIRXS	CLSWAIRXS
SM	4.5 - 7	6 - 8	CLTWAIRSM	CLSWAIRSM
MD	7.5 - 10	8.5 - 11.5	CLTWAIRMD	CLSWAIRMD
LG	10 - 13	11 - 15	CLTWAIRLG	CLSWAIRLG
XL	13+	15+	CLTWAIRXL	CLSWAIRXL



CORELINE NON-PNEUMATIC WALKER W/ METAL UPRIGHTS

L-Code: L4387 / L4386

- Provides maximum comfort and stability for fractures, soft tissue injuries or severe sprain of the foot and ankle. Metal uprights can be bent to accommodate ankle edema or bulky dressings.
- Low natural heel height, with rocker bottom and slight toe flexion promotes natural gait

Size	Men's Shoe Size	Women's Shoe Size	Tall Product No.	Short Product No.
XS	2 - 4	3 - 5	CL0020	CL0010
SM	4.5 - 6	5.5 - 7	CL0022	CL0012
MD	6.5 - 10	7.5 - 11	CL0024	CL0014
LG	10.5 - 12	11.5 - 13	CL0026	CL0016
XL	12 - 17	13+	CL0028	CL0018



CORELINE PNEUMATIC WALKER W/ METAL UPRIGHTS

L-Code: L4361 / L4360

- Has easily adjustable air cells to increase stability and compression while decreasing pain and swelling
- High and low options available

Size	Men's Shoe Size	Women's Shoe Size	Tall Product No.	Short Product No.
XS	2 - 4	3 - 5	CL0040	CL0030
SM	4.5 - 6	5.5 - 7	CL0042	CL0032
MD	6.5 - 10	7.5 - 11	CL0044	CL0034
LG	10.5 - 12	11.5 - 13	CL0046	CL0036
XL	12 - 17	13+	CL0048	CL0038

NEW



CORELINE MULTI PODUS BOOT

L-Code: L4397 / L4396

- Off-loading heel to help prevent ulcers
- Holds foot at a neutral position to help prevent contractures
- Adjustable kickstand
- EVA sole for ambulation

Size	with Fleece	without Fleece
REGULAR	CLMPBFL-REG	CLMPB-REG
XL	CLMPBFL-XL	CLMPB-XL

Knee Bracing



CORELINE ROM WRAP KNEE UNIVERSAL

L-Code: L1833 / L1832

- Malleable aluminum hinge allowing for a more custom fit.
- Dual axis motion control hinges provide dial in adjustments with flexion/extension stops at 0, 15, 30, 60 and 90 degrees.
- Front opening for easy on/off and custom

Size	Product No.
12"	CLROMWAKUNV-12
16"	CLROMWAKUNV-16



CORELINE WRAP AROUND HINGED KNEE

L-Code: L1820

- Neoprene material
- Open POP with condyle pad

Size	Measurement	w/ Condyle Pad	w/o Condyle Pad*
XS	13.5" - 15.5"	CL0152	CL0162
SM	15.5" - 18"	CL0153	CL0163
MD	18" - 20.5"	CL0154	CL0164
LG	20.5" - 23"	CL0155	CL0165
XL	23" - 25.5"	CL0156	CL0166
2XL	25.5" - 28"	CL0157	CL0167
3XL	28" - 30.5"	CL0158	CL0168



CORELINE OA UNLOADER

L-Code: L1851 / L1843

- Low profile, lightweight design
- Single upright
- Easy to apply push buttons

Size	Product No.
UNIVERSAL LEFT	CLOAKBUNV-L
UNIVERSAL RIGHT	CLOAKBUNV-R




CORELINE TRI PANEL KNEE IMMOBILIZER

L-Code: L1830

- Adjustable three piece design accommodates multiple circumferences
- Lightweight aluminum posterior, medial and lateral stays are malleable to conform to the contours of the knee

Size	w/o Condyle Pad*
12"	CLKI12IN
16"	CLKI16IN
20"	CLKI20IN
22"	CLKI22IN
24"	CLKI24IN



CORELINE UNIVERSAL POST OP BRACE

L-Code: L1833 / L1832

- Lightweight, low profile design
- ROM hinge with drop lock
- Fully customizable
- One size fits all

Size	Product No.
UNIVERSAL	CLPOKBUNV



CORELINE POSTERIOR NIGHT SPLINT

L-Code: L4397 / L4396

- Adjustable dorsiflexion straps provide a gentle stretching of the plantar fascia and Achilles tendon
- Soft lined foam for increased comfort
- Slip-resistant tread for increased user safety

Size	Men's Shoe Size	Women's Shoe size	Product No.
SM	<7	<7.5	CLPNSSM
MD	7.5 - 9	8 - 9.5	CLPNSMD
LG	9.5 - 11.5	10 - 12.5	CLPNSLG



CORELINE ANTERIOR NIGHT SPLINT

L-Code: L4397 / L4396

- Adjustable dorsiflexion straps provide a gentle stretching of the plantar fascia and Achilles tendon
- Soft lined foam for increased comfort
- Easy to use, low profile design

Size	Men's Shoe Size	Women's Shoe size	Product No.
SM / MD	<9	<9.5	CLANSSM/MD
LG / XL	9.5<	10<	CLANSLG/XL



CORELINE FIGURE 8 ANKLE

L-Code: L1902

- Nylon stabilizing straps restrict inversion or eversion ankle motion
- Fits either the left or the right foot

Size	Ankle Circumference	Product No.
XS	10" - 11"	CLF8AXS
SM	11" - 12"	CLF8ASM
MD	12" - 13"	CLF8AMD
LG	13" - 14"	CLF8ALG
XL	15" - 16"	CLF8AXL



CORELINE AIR/GEL STIRRUP ANKLE BRACE

L-Code: L4350

- Designed to protect a sprained ankle while enabling healing and continued use by the patient
- Ankle wrap ensures control of swelling to speed recovery while minimizing ankle rotation
- Air/Gel bladder may be used for cold therapy to treat edema

Size	Height	Product No.
XS	6"	CL0090
SM	8.75"	CL0092
MD	9.5"	CL0094
LG	10"	CL0096



CORELINE THUMB SPICA

L-Code: L3809 / L3807

- Durable neoprene construction has removable aluminum stays to support the palmar surface of wrist and extensor surface of the thumb
- Provides wrist immobilization while allowing full hand functions

Size	Wrist Circ.	Left Product No.	Right Product No.
SM	4.75" - 5.5"	CLTSLSM	CLTSRSM
MD	5.5" - 6.25"	CLTSLMD	CLTSRMD
LG	6.25" - 7"	CLTSLLG	CLTSRLG
XL	7" - 7.75"	CLTSLXL	CLTSRXL
UNIVERSAL LEFT		CLTSUNIVL	
UNIVERSAL RIGHT		CLTSUNIVR	



CORELINE COCK UP WRIST SPLINT

L-Code: L3908

- Non-removable anatomically designed aluminum dorsal stay maintains positioning with padding for comfort
- Designed to fit below palmar crease for full finger dexterity

Size	Wrist Circ.	Left Product No.	Right Product No.
SM	4.75" - 5.5"	CLWSLSM	CLWSRSM
MD	5.5" - 6.25"	CLWSLMD	CLWSRMD
LG	6.25" - 7"	CLWSLLG	CLWSRLG
XL	7" - 7.75"	CLWSLXL	CLWSRXL
LEFT UNIVERSAL		CLWSUNIVL	
RIGHT UNIVERSAL		CLWSUNIVR	



CORELINE CERVICAL COLLAR / 3" AND 4"

L-Code: L0120

- Firm density foam collar designed to support cervical spine in neutral position
- Available in 3" or 4" to fit a wide range of patient neck sizes

Size	Neck Cir.	3" Product No.	4" Product No.
XS	9"-13"	CL0100	—
SM	11"-16"	CL0102	CL0112
MD	13"-18"	CL0104	CL0114
LG	15"-20"	CL0106	CL0116
XL	16"-22"	CL0108	—
UNIV	14"-25"	CL0120	CL0125



CORELINE SHOULDER PAD UNIVERSAL

L-Code: L3670

- Sling has multiple adjustments to accommodate all body types
- Large pillow provides abduction with comfort
- Waist strap that provides extra stability

Size	Product No.
UNIVERSAL	CLSPUNV



CORELINE L0627 BACK BRACE SIZE UNIVERSAL

L-Code: L0642 / L0627

- One Size fits all, low profile design
- Pulley System comfortably tightens brace for maximum comfort

Size

UNIVERSAL

Product No.

CLLSO100



CORELINE L0631 BACK BRACE SIZE UNIVERSAL

L-Code: L0648 / L0631

- One Size fits all, low profile design
- Pulley System comfortably tightens brace for maximum comfort

Size

UNIVERSAL

Product No.

CLLSO200



CORELINE L0637 BACK BRACE SIZE UNIVERSAL

L-Code: L0650 / L0637

- One Size fits all, low profile design
- Pulley System comfortably tightens brace for maximum comfort
- Removable side panels that provide maximum support

Size

UNIVERSAL

Product No.

CLLSO300

NEW



CORELINE EXTENSION PANEL

- Sold as pair
- Measures 12" in length

Size

UNIVERSAL

Product No.

CLEXTTP

NEW



CORELINE L0456/L0457 UNIVERSAL TLSO

L-Code: L0456 / L0457

- Effectively assist postural correction and re-education
- Provides rotational stability and motion restriction
- Adjustable chest strap for a custom fit
- PDAC approved
- Low profile design featuring breathable mesh belt and moisture-wicking lining
- No grip closure for patients lacking hand strength
- Patient friendly design

Size

UNIVERSAL (SX - 3XL); (22" - 55")

Product No.

CLTLSO56

NEW



CORELINE L0464 UNIVERSAL TLSO

L-Code: L0464

- Provides optimal motion restriction
- Removable sternal pad
- PDAC approved
- Low profile design featuring breathable mesh belt and moisture-wicking lining
- No grip closure for patients lacking hand strength
- Patient friendly design

Size

UNIVERSAL (SX - 3XL); (22" - 55")

Product No.


CLTLSO64



CORELINE CRUTCHES

- 300 lb Weight Limit

Size	Product No.
Youth 4' 6" - 5' 2"	CLCR-YTH
Adult Regular 5' 2" - 5' 10"	CLCR-REG
Adult Tall 5' 10" - 6' 6"	CLCR-TALL





800.288.5570
PUNCHOUT TO HANGER STORE

April 2023